



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

## ZAP™ 1200 uL Aerosol Filter Pipet Tips for Rainin® LTS Pipettors, in 96 Racks, Sterile

<b>Fluid Measure</b>	μL	<b>Freezable</b>	NA
<b>Fluid Volume</b>	1200	<b>Freeze Temp</b>	
<b>Height (mm)</b>	88.9	<b>Autoclavable</b>	121°C 15min
<b>Width (mm)</b>		<b>Sterile</b>	Yes
<b>Length (mm)</b>		<b>Sterile Processing</b>	Certified as sterile with a sterility assurance level (SAL) of 10 <sup>-6</sup>
<b>Pack Qty</b>	576	<b>Centrifugation</b>	NA
<b>Case Qty</b>	4608	<b>Centrifugation Speed (RCF)</b>	
<b>Package Type</b>	96 Rack	<b>Shelf Life (Months)</b>	36
<b>Material</b>	Polypropylene	<b>Certification</b>	
<b>Cap Material</b>		<b>Regulatory</b>	

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

**Additional Testing:** Protease, Human DNA, ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)



This Technical Data Sheet is provided as a service and for documentation only. Labcon North America retains all rights to the information, procedures, and results of the material contained hereon as proprietary information. Any copying and distribution of this technical data without prior written authorization is prohibited by law. The use of this information and documentation does not discharge the user from their responsibility to test and use this product appropriately. This item is for use solely by laboratory professionals. This document is automatically generated and is © **Copyright 10-May-2021** . All Rights Reserved. Technical Data Sheet Version 1/2020