# GM70 Hand-Held Carbon Dioxide Meter for Spot-Checking Applications



The Vaisala CARBOCAP\* Hand-Held Carbon Dioxide Meter GM70 is the demanding professional's choice for hand-held carbon dioxide measurement. The meter consists of the indicator (center) and probe, used either with the handle (left) or pump (right).

The Vaisala CARBOCAP® Hand-Held Carbon Dioxide Meter GM70 is a user-friendly meter for demanding spot measurements in laboratories, greenhouses and mushroom farms. The meter can also be used in HVAC and industrial applications, and as a tool for checking fixed CO<sub>2</sub> instruments.

# Vaisala CARBOCAP® technology


The GM70 incorporates the advanded Vaisala CARBOCAP\* Sensor. The patented sensor has unique reference measurement capabilities. Its critical parts are made of silicon; this gives the sensor outstanding stability over both time and temperature. In addition, the measurement accuracy of the sensor is not affected by dust, water vapor or most chemicals.

The GM70 has a two-year recommended calibration interval.

# Two sampling methods

The handle is for hand-held diffusion sampling. The GM70 pump enables pump-aspirated sampling from locations difficult to access otherwise. It is also ideal for comparisons with fixed  $\mathrm{CO}_2$  transmitters.

Due to its short warm-up time, the GM70 is ready for use almost immediately. The menu-based interface is easy to use, and a graphical LCD display with data logger function makes the handling of measurement data convenient. The optional MI70 Link Windows\* software in



GM70 displays  $CO_2$  concentration both numerically and graphically.

## **Features/Benefits**

- Proven Vaisala CARBOCAP\* reliability
- Two optional sampling methods: diffusion or pump aspiration
- User-friendly meter with multilingual user interface
- Numerical and graphical display of measurements
- Data can be logged and transferred to PC via MI70 Link software
- Wide selection of measurement ranges
- Easy recalibration using the interchangeable probes
- Suitable for field checking of fixed CO<sub>2</sub> instruments
- Short warm-up time
- Compact and versatile

combination with a USB connection cable provides an easy way to handle data in a PC environment.

#### **Interchangeable probes**

The GM70 uses the same probes as Vaisala CARBOCAP° Carbon Dioxide Transmitter Series GMT220 fixed  $\mathrm{CO}_2$  instruments. By plugging different probes into the handle or pump, the user can easily change the measurement range.

#### Two inputs provide flexibility

The GM70 has two probe inputs. This allows differences to be measured, and is also a practical feature when checking probes against a calibrated probe. Vaisala's relative humidity and dewpoint probes can also be used simultaneously.

The meter can also be used as a calibration check instrument for Vaisala's GM20 and GM220 series fixed  ${\rm CO}_2$  instruments. GMP220 probes can even be adjusted by using the GM70 meter.

GM70 CARBON DIOXIDE

# **Technical Data**

#### CO. volume concentration measurement

co, volume concentration measurement		
Measurement ranges		
High concentrations	0 2 %	
short probe (GMP221)	0 3 %	
•	0 5 %	
	0 10 %	
	0 20 %	
Low concentrations	0 2000 ppm	
long probe (GMP222)	0 3000 ppm	
	0 5000 ppm	
	0 7000 ppm	
	0 10,000 ppm	
Accuracy (including repeatability, n	on-linearity and calibration	
uncertainty) at 25 °C and 1013 hPa	·	
GMP221	$\pm (1.5\% \text{ of range} + 2\% \text{ of reading})$	
(applies for concentrations above 29	% of full scale)	
GMP222	$\pm (1.5\% \text{ of range} + 2\% \text{ of reading})$	
Temperature dependence, typical	-0.3 % of reading / °C	
Pressure dependence, typical	+0.15% of reading/hPa	
Long-term stability	<±5 %FS/2 years	
Response time (63 %)	•	
GMP221	20 seconds	
GMP222	30 seconds	

# **Measurement environment**

Temperature  $\begin{array}{ccc} -20 \dots +60 \text{ °C } (-4 \dots +140 \text{ °F}) \\ \text{Relative humidity} & 0 \dots 100 \text{ %RH non-condensing} \\ \text{Operation pressure} & 700 \dots 1300 \text{ hPa} \\ \text{Flow range (diffusion sampling)} & 0 \dots 10 \text{ m/s} \\ \end{array}$ 

30 seconds, 15 minute full specifications

#### Probe, handle & pump general

riobe, namale & painp g	Cilciai
Sensor	Vaisala CARBOCAP®
Housing material	
GMP221/222 probe	PC plastic
GMH70 handle	ABS/PC blend
GM70 Pump	aluminium IP54
Storage temperature	-30 +70 °C (-22 +158 °F)
Storage humidity	0 100 %RH non-condensing
Weight	0
GMH70 with GMP221/222 probe	230 g
GM70 Pump with GMP221/222 prob	oe 700 g
*	

## **Dimensions**

Warm-up time

GMH70 handle with probe	
length, short probe (GMP221)	270 mm (10.6 in.)
length, long probe (GMP222)	350 mm (13.8 in.)
Space required for GM70 pump	
length (with long probe & cable)	300 mm (12 in.)
width	85 mm (3.3 in.)
height	36 mm (1.4 in.)
MI70 indicator	
length	200 mm (7.87 in.)
width	80 mm (3.1 in.)
height	43 mm (1.7 in.)



35 Vantage Point Drive Rochester, NY 14624 Call 1.800.800.5001

# MI70 Indicator general

IVIII/ O IIIIGICO	ator gener	мі
Menu languages	English, (	Chinese, French, Spanish, German,
	J:	apanese, Russian, Swedish, Finnish
Display	LCD with	backlight, graphic trend display of
	any param	eter, character height up to 16 mm
Probe inputs		1 or 2
Power supply	Rechargeable Ni	MH battery pack with AC-POWER
Analog output	_	0 1 VDC
Output resolution	1	0.6 mV
Data interface	RS232 (access	sible only with MI70 Link software)
Data logging capa	acity	2700 points
Alarm	•	Audible alarm function
Operating temper	rature range	-10 +40 °C (+14 +104 °F)
Operating humid	ity range	non-condensing
Housing material		ABS/PC blend
Housing classifica	ation	IP54
Weight		400 g

# **Battery operation time**

Continuous use	
with handle	better than 8h at +20 °C (+68 °F)
with pump	better than 5h at +20 °C (+68 °F) without load
Data logging use	up to 30 days

# **Electromagnetic compatibility**

EN 61326-1, Portable Equipment.

#### **Accessories**

Connection cable for fixed CO <sub>2</sub> instruments		
GMT220, GMM220, GMD20 and GMW20	GMA70	
USB cable for Mi70	219687	
Analog output cable for 0 1 VDC analog output	27168ZZ	
Windows® software, includes PC connection cable	e MI70LINK	
Calibration adapter	26150GM	
Carrying case	MI70CASE	
Battery, NiMH 4.8V	26755	
Spare probe	GMP221, GMP222	
(use the order form to define measurement range etc.)		

Nafion Membrane Tubing 212807GM



The GM70 in pump-aspirated mode is used for example to check incubators in the field.

 $CARBOCAP^{\ast}$  is a registered trademark of Vaisala. Specifications subject to change without prior notice. ©Vaisala Oyj

